

## **MICROBIOLOGY STAINING PROCEDURES – QUICK NOTES**

### **1. Simple Staining**

- Purpose: To observe shape, size, arrangement of bacteria.
- Dyes used: Methylene blue, Crystal violet, Safranin.
- Steps:
  1. Prepare smear → air dry → heat fix.
  2. Add basic dye for 1 minute.
  3. Rinse with water & blot dry.

### **2. Gram Staining (Differential Stain)**

- Purpose: To differentiate Gram-positive & Gram-negative bacteria.
- Principle: Difference in cell wall peptidoglycan thickness.
- Steps:
  1. Primary stain: Crystal violet (1 min).
  2. Mordant: Iodine (1 min).
  3. Decolorizer: Alcohol/Acetone (few seconds).
  4. Counterstain: Safranin (1 min).
- Results:

Gram + = Purple  
Gram – = Pink

### **3. Acid-Fast Staining (Ziehl–Neelsen)**

- Purpose: To identify Mycobacteria (waxy cell wall with mycolic acid).
- Steps:
  1. Primary stain: Carbol fuchsin with heat.
  2. Decolorizer: Acid alcohol.
  3. Counterstain: Methylene blue.
- Results:

Acid-fast = Red  
Non–acid-fast = Blue

### **4. Endospore Staining (Schaeffer–Fulton)**

- Purpose: To detect bacterial spores (e.g., Bacillus, Clostridium).
- Steps:
  1. Primary stain: Malachite green with heat.
  2. Rinse with water.
  3. Counterstain: Safranin.
- Results:

Spores = Green  
Vegetative cells = Red

### **5. Capsule Staining (Negative Staining)**

- Purpose: To visualize bacterial capsules (non-ionic & non-staining).
- Steps:
  1. Mix organism in India ink / Nigrosin.
  2. Air dry (no heat fixing).
  3. Counterstain: Crystal violet.
- Results:

Capsule = Clear halo  
Background = Dark  
Cells = Purple

### **6. Flagella Staining (Hanging Drop)**

- Purpose: To observe motility.
- Principle: Use of mordant to thicken flagella for visibility.

**\*\*\*End of Notes\*\*\***